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INSTRUCTIONS FOR USE

A) PRODUCT IDENTIFICATION:

TECHNICAL NAME: Surgical instruments for dental implants.

B) SOCIAL REASON AND ADDRESS OF THE MANUFACTURER:

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GINO KOPP

CRO No. 7679- PR

C) INDICATION, PURPOSE, USE, APPLICATION AND PRINCIPLE OF PRODUCTS:

The instruments are used and indicated for dental surgery.

“NON-STERILE PRODUCT - REUSABLE”

D) INFORMATION ON USE:

D-1) HOW TO USE:

After the correct sterilization of the components, they must be handled aseptically and according to standards of biosafety.

ADVERSE EFFECTS:

The instruments are used and indicated for dental surgery, that is, adverse effects only will occur if the choice of the instrument is improper, and may cause damage to the component to be installed.

D-2) ADVERSE REACTIONS:

No adverse reactions to the product have been reported, as the components of the instruments are made from compatible materials.

E) PRECAUTIONS, CARES AND CLARIFICATIONS ABOUT POSSIBLE RISKS:

E-1) CONDITIONS OF USE:

NON-STERILE INSTRUMENTALS MUST BE PROPERLY STERILIZED BEFORE USE.

RESTRICTED USE ON QUALIFIED PROFESSIONALS.

It is necessary to take special care and its use should only be done by a professional trained to perform cleaning, asepsis, descaling, washing and sterilization procedures for instruments subject to contamination and that are not disposable. Use according to sterilization procedures.

All instruments must be discarded and no longer re-sterilized when they have scratches, cracks or dents high intensity or when it loses its cutting action.

E-2) GENERAL PRECAUTIONS:

The instruments require specialized surgical procedures, should only be used by dental surgeons enabled, including: diagnosis, preoperative planning and surgical protocol. Using the product without knowledge appropriate techniques and / or procedures and conditions, may harm the patient leading to unsatisfactory results.

E-3) STORAGE AND TRANSPORT:

Keep the product at room temperature (15°C to 30 ° C), away from light and heat. The instruments have a case that helps to protect against falls and collisions, however all must be transported in a safe, clean and closed vehicle.

F) INFORMATION ON REUSE PROCEDURES:

The non-sterile instruments that make up the Kit are subject to sterilization and reuse and must be strictly followed the washing or descaling and sterilization procedures, as described below:

Washing or descaling: removing organic matter from the instruments. Must be started more quickly possible after use in dental implant installation surgeries.

a) Always use deionized water for this procedure;

If the water is heated, to facilitate this cleaning step, the temperature must be between 30°C and 40°C;

Use 10% ENZYMATIC detergent, such as Prolystica. Dilution should be carried out as recommended manufacturer;

b) Immerse the entire component body in the detergent solution, keeping the solution in contact with the component for at least 3 minutes, if the chosen detergent requires more time, this must be respected;

c) Rub each component with a soft bristle brush, at least 5 times, from the proximal to the distal direction. Repeat this procedure until the visible dirt is eliminated, making sure that all the recesses have been washed;

d) Rub the inner surface of each lumen with a soft brush, adjusted to the size of the lumen, at least 5 from the proximal to the distal direction. Repeat this procedure until the visible dirt is removed;

e) Never use steel wool or sponges and / or abrasive products, as this may damage the instruments in use;

f) Do not accumulate instruments in large quantities, on top of each other, to prevent deformation of parts smaller and delicate and thus also do not scratch the polished surfaces. Always handle a few parts at a time;

g) Rinse the component in deionized water for at least 1 minute;

h) Dry the instrument with a soft and clean cloth.

IMPORTANT! Always use deionized water to wash the stainless steel. Tap water contains chlorine and causes oxidation of stainless steel. Do not use descaling products, as they tend to darken the parts and cause oxidation. The enzymatic liquid in a concentration greater than 10% favors oxidation, as well as the drying of parts containing residues of the solution of cleaning. The drying of the pieces is extremely important before storage and sterilization, because the accumulation of moisture in the products is harmful and can cause oxidation spots.

Sterilization: it is the procedure that aims at the total elimination of microorganisms. Use sterilization procedure by **SATURATED STEAM UNDER AUTOCLAVE PRESSURE.**

- a) Wrap individually in packages, according to the sterilization process. Preferably use a film envelope and surgical grade paper, suitable for steam sterilization;
- b) Do not open the autoclave prematurely, to avoid rapid condensation;
- c) Do not open the autoclave quickly, letting all the steam out, before the drying cycle is complete;
- d) In an autoclave the instrument must remain for 30 minutes at a temperature of 121°C;
- e) The drying time is 15 minutes;
- f) Clean the autoclave rigorously and periodically, removing dirt and eventual excess rust.

Important:

Never store the material without complete drying. Sterilize the day before or the day of the procedure.

Note:

During the sterilization process, no temperature controller may have a reading above 126° C in the first 60 seconds and after this time the minimum allowed variation must equal the sterilization temperature 121° C and the maximum of 124° C.

These washing and sterilization procedures comply with the requirements specified in the standards: **ISO 11737-1; ISO 11737-2; ISO 17665-1; DRC 16: 2013; EN ISO 13485: 2016; ISO 17664: 2017.**

I) DISPOSAL OF MATERIALS:

All consumables used in surgery to install dental implants and their components should be disposed of in hospital waste.

J) BATCH, MANUFACTURE AND VALIDITY:

The product has an indefinite expiration date.

H) REGISTRATION WITH ANVISA Nº:

Product	Code ANVISA
KOPP TACK KIT.....	80204110002
KOPP PROSTHETIC KIT	80204110004
KOPP DISPOSABLE INSTRUMENTAL KIT	80204110005
KOPP ELEVATOR KIT	80204110006
KOPP SURGICAL KIT	80204110013
KOPP MANUAL BONE GRINDERS	80204110014
KOPP EXPANSOR KIT	80204110015
ORTHODONTIC KIT / FIXATION	80204110016
KOPP IMPLANT MINI KIT	80204110019
KOPP GUIDE KIT	80204110020
KOPP FIXING KIT	80204110021
KOPP PROSTHETIC SELECTION KIT	80204110022
MINI SURGICAL KIT KOPP	80204110023
KOPP GUIDE HIN KIT	80204110024

Curitiba, November 11, 2020

LEGAL REPRESENTATIVE / TECHNICAL RESPONSIBLE

GINO KOPP
CPF: 552.626.519-68
CRO Nº 7679 - PR