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INDICATIONS FOR USE:

Kopp CM Screw Implant is a dental implant devised to restore surgically a patient's chewing function. It can be placed in the bone of the upper or lower jaw to provide a substructure support for the prosthetic devices, therefore providing an artificial root means for a single or multiple unit prosthetic appliance. Dental implants can be placed in a conventional two stage surgical process with an option for transmucosal healing or a single stage surgical process for immediate loading. In this immediate clinical case, the primary stability needs to be in a good standpoint and with an appropriate occlusal loading. Kopp CM Screw can be placed in any bone density.

INFORMATION ABOUT THE USE:

In a preoperative evaluation, the patient treatment must be evaluated in order to determine the appropriate specifications of the products to be used in the surgery. Complimentary, the most suitable technique will be chosen by the dentist as well, according to the dental surgeon knowledge. It is highly recommended to prepare a controlled protocol of the postoperative period, including actions taken by the surgeon in a case of unexpected intervention during the surgical procedure. Another important documentation to be attached on the patient case is the clinical follow-up evaluation report, containing the surgical results observed by the clinician. Furthermore, the prophylactic treatment should be adopted in the preoperative period as well as in the immediate and late postoperative period as stated in the classical literature.

Dental implants surgical interventions should not be considered in circumstances when patients are not healthily apt in a medical standpoint. A preoperative clinical evaluation should be carefully conducted by the oral surgeon, especially on patients who demonstrates specific concerns or systematic factors that might interfere in the healing process of the bone or soft tissues. Another important factor to be analyzed beforehand and during the surgery is the appropriate volume of bone in the surgical site where the dental implant will be placed. In order to achieve a satisfactory treatment, the patient must have an appropriate quantity of residual bone in the surgical site and the oral surgeon must select the correct implant dimension. This clinical evaluation is extremely important for the long-term preservation of the product functionality. It prognosticates the functionals mechanism to support the patients' occlusion. Therefore, proportioning a satisfactory result after the surgical treatment if the adequate prognostic and technique are applied by the Implantologist. In cases of not favorable bone volume, enlargement bone technics must be considered. Otherwise, in any other circumstances where the professional is uncertain about the adequate bone volume, the procedure should not be executed.

*In any case where the oral surgeon requires further analyses of a specific product, the adhesive label of the product must be attached to the patient's records documents, aiming at a satisfactory case analyzes and adequate traceability. Otherwise, the company will be exonerated of any legal responsibility towards any concerns regarding to the product.

All of the documentation mentioned above should be provided such as the preoperative planning, the follow-up evaluation report, containing the surgical complications and the risk that might occur during the dental implant surgery. On the case of installed prosthesis, the implant supported prosthesis should be provided. In addition, another requirement is the adhesive label, which discriminates the product batch, material and date of manufactured and name of the product as well. In addition, a radiographic follow-up must be provided along with the patients' health state condition of life. Furthermore, the surgeon must consider that there is a low-probability percentage of osseointegration failure due to various factors stated in the dentistry literature. Therefore, in any of these occasions, the company will not take the responsibility. Another important point that should be stated by the oral surgeon, it is that, dental implant and prosthesis are not ensured to last for a long period of time due to various aspect and are scientific proven. Thus, no warranties should be given regarding those circumstances. Patients must be warned about potential paresthesia, dysesthesias and risk of infections. Moreover, low bone quality and quantity, poor hygiene, infections, systematic disturbances can cause the loss of the implant in the future. In case of any damage of the packaging or any damage caused in the product, the dental implant must be discarded immediately and not utilized in any circumstances in the procedure. According to Professor Bränemark, in order to maintain and obtain a successful osteointegration, the oral surgical products must be utilized by professionals duly qualified to apply the adequate method of application.

IMPLANT MEASURES COMMERCIALY AVAILABLE (IMPLANT WIDTH AND DIAMETER)

Implant CM Screw Platinum

| | 7.0 mm | 8.0 mm | 9.0 mm | 11.0 mm | 13.0 mm | 15.0 mm |
|-----------|---------|---------|---------|---------|---------|---------|
| ø 3.5 mm | 140.714 | 140.715 | 140.716 | 140.717 | 140.718 | 140.719 |
| ø 3.75 mm | 140.720 | 140.721 | 140.722 | 140.723 | 140.724 | 140.725 |
| ø 4.0 mm | 140.726 | 140.727 | 140.728 | 140.729 | 140.730 | 140.731 |
| ø 4.3 mm | 140.732 | 140.733 | 140.734 | 140.735 | 140.736 | 140.737 |
| ø 5.0 mm | 140.738 | 140.739 | 140.740 | 140.741 | 140.742 | 140.743 |

Implant CM Screw Standard

| | 7.0 mm | 8.0 mm | 9.0 mm | 11.0 mm | 13.0 mm | 15.0 mm |
|-----------|---------|---------|---------|---------|---------|---------|
| ø 3.5 mm | 140.744 | 140.745 | 140.746 | 140.747 | 140.748 | 140.749 |
| ø 3.75 mm | 140.750 | 140.751 | 140.752 | 140.753 | 140.754 | 140.755 |
| ø 4.0 mm | 140.756 | 140.757 | 140.758 | 140.759 | 140.760 | 140.761 |
| ø 4.3 mm | 140.762 | 140.763 | 140.764 | 140.765 | 140.766 | 140.767 |
| ø 5.0 mm | 140.768 | 140.769 | - | - | - | - |

Implant CM Screw Slim

| | 7.0 mm | 8.0 mm | 9.0 mm | 11.0 mm | 13.0 mm | 15.0 mm |
|----------|--------|--------|---------|---------|---------|---------|
| ø 3.0 mm | - | - | 140.139 | 140.140 | 140.141 | - |

WARNING:

Only qualified implantology professionals who have dental implant surgical knowledge, regarding techniques and procedures, are allowed to utilize the products, adhering to the appropriate application and utilization of the product. In addition, it is extremely important to respect legal requisites to be able to practice the surgery. Otherwise, the use of this product is prohibited. Furthermore, the lack of dental implant surgical knowledge can cause harm to patient, leading to non-satisfactory results.

HOW TO USE:

- After the diagnosis, a surgical planning must be done, determining where the implant will be placed and adequate dental implant measures to be adopted.
- Administer the preoperative medication.
- Anesthetize the patient and subsequently make the incision.
- Utilize the necessary sequence of drills depending on the implant diameter selected by the oral surgeon.

RECOMMENDED DRILLING SEQUENCE

| | Ø 3.0 | Ø 3.5 | Ø 3.75 | Ø 4.0 | Ø 4.3 | Ø 5.0 |
|-------|-------|-------|--------|-------|-------|-------|
| lança | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2.0 | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2.8 | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 3.0 | | ✓ | ✓ | ✓ | ✓ | ✓ |
| 3.3 | | | ✓ | ✓ | ✓ | ✓ |
| 3.8 | | | | | ✓ | ✓ |
| 4.3 | | | | | | ✓ |

- Do not utilize drilling rotation higher than 900 RPM and always irrigate the surgical area with physiological serum.
- Be extremely cautious in anatomical sites that demonstrate vital vessels or nerves, respecting the sensitive and nutritious areas where the implant will be placed. Do not aspirate the clot.
- After drilling the surgical alveolus, open the primary and secondary packagings.
- Remove the dental implant from inside of the cocoon with the Monta key instrumental docked by the contra-angel. In order to hold the dental implant safely, it is important to import a pressure with the key towards the internal hexagon of the dental implant. It is important to not over pressure the force when doing it, otherwise it can cause the damage of the product.
- Take out the dental implant cautiously from the cocoon in a 180-degree position. It is also important to check if there is retention in the key.
- Carry the implant in the surgical site to insert it.
- Screw drive the dental implant with the Torque meter Ratchet until the instrumental indicates 45 Newtons. It is extremely important to not exceed this parameter.
- Counter angle up until the point where the platform reaches the subosseous position, proportionating an appropriate depth according to the selected implant.

- Place the implant cap or healing component, selecting the appropriate one conforming the implant measures.
- Suture and dispose instructions to the patient regarding the postoperative care.
- Follow up the process of healing and osseointegration.

POSITION FOR SETTING THE DENTAL IMPLANT:

The implant platform must be set 2mm below the osseous level in order to allow the most suitable emergency profile and to respect the biological surgical site. The implant of the maximum diameter and length, as allowed by the structures of the bones, must be chosen by the oral surgeon. The use of this product is not recommended for patients who demonstrates the following clinical statuses:

- Non-controlled disturbances of the metabolism (for instance: diabetes mellitus, hyperparathyroidism).
 - Generalized renal dysfunction and serious liver diseases.
 - Presence of vascular damages or impairment at the area of the implant.
 - Osteomyelitis at the area of the surgery.
 - Patients who received irradiation a short time ago.
 - Patients with sensitiveness to the titanium.
 - Patients with generalized infection.
 - Pregnant and lactating women (as per medical criterion).
 - Smokers, drunkards and chemical dependents (as per professional criterion).
 - Patients with poor mouth hygiene.
 - In cases of excessive occlusal function.
 - Set the torque 60N.
 - In case that the locking of the Implant Assembler occurs during the setting, turn the ratchet in the opposite direction by a gentle counter torque.
- No adverse reactions to the product were reported as the dental implants, made of titanium, are biocompatible. Considering that the product is commercialized in the sterile, the professional will be responsible for handling it by following asepsis and biosafety criteria. Material should not to be sterilized at the consultation room.

ADVERSE EFFECTS:

Precautions, Cautions and Clarifications about Possible Risks: Sterile; Single Use Only; Used Restricted to Qualified Professionals; Sterilized by Using Ethylene Oxide.
Do not Utilize if There is any Damage Signal on the Packaging.

GENERAL PRECAUTIONS:

The surgical products and the prosthetic restoration products utilized to get and to keep the osseointegration, as described by Professor Brånemark et al, must be utilized by a professional duly qualified for the respective method of application. Each patient will have to be carefully examined and evaluated, viewing the determination of the radiographic, psychological and physical status as well as the dental deficits regarding the bone or the adjacent soft tissue that can influence upon the final result of the surgical intervention.

STORAGE AND TRANSPORTATION:

The product must be kept in the ambient temperature, away from light and heat and should not be utilized after the validate expiration date or if the package is damaged. Please, keep away from children.

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