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INDICATIONS FOR USE:

The Cone Morse Screw is a dental implant developed to surgically reestablish a patient's masticatory function. Kopp's intermediates (abutment, component) on implants enable prosthetic rehabilitation safely. At first, dental implants can be placed in a conventional two-stage surgical procedure, or a single stage can be chosen by the dentist. The Cone Morse Screw can be placed in any type of bone density, cortical or medullary, upper or lower jaw.

Note: In cases of immediate loading, it is necessary for the dentist to prove that there is the primary stability and torque required in the installation of the dental implant for the installation of the cosmetic prosthesis, without occlusal loading.

DESCRIPTION:

Cone Morse Screw dental implants are produced in pure titanium (grade IV) and available in various heights and diameters on the market. All Kopp dental implants are marketed with the Surface Speed - rough surface, obtained by double treatment: abrasive blasting and acid subtraction. Its macrogeometry stands out due to the following characteristics:

- Morse cone prosthetic interface with internal hexagonal indexer;
- Biological chamfer and conical middle region;
- Double threads of trapezoidal profile;
- Active tip that favors its high penetration;
- Ability to understand bone during installation.

CONTRAINDICATIONS:

The contraindications to the use of dental implants are related to the appropriate selection of the patient according to general dental and medical restrictions. These include, but are not limited to:

- Coagulation disorders.
- Disturbance of healing or bone regeneration.
- Occlusal disorders, as well as small interocclusal spaces.
- Poor oral hygiene.
- Infection or inflammation in the oral cavity.
- Patient with insufficient bone volume or soft tissue coverage.
- Patient who is not prepared to undergo total dental arch rehabilitation.
- Untreated screws.
- Immune suppression therapy.

SIDE EFFECTS:

Possible side effects of installing dental implants can be:

- Temporary anesthesia or paresthesia and masticatory impediment.
- Temporary local edema, edema and hematoma.

MAGNETIC RESONANCE IMAGE

The safety and compatibility of Kopp dental implants with the MRI environment have not been evaluated. The heating, distortion or displacement suffered by Kopp dental implants in the MRI environment has not been tested. The safety of Kopp dental implants in the MRI environment is unknown. Performing an MRI on a patient with this device can result in damage to the patient.

COMMERCIALLY AVAILABLE IMPLANT MEASURES (IMPLANT LENGTH AND DIAMETER)

Implant CM Screw Platinum

	7.0 mm	8.0 mm	9.0 mm	11.0 mm	13.0 mm	15.0 mm
ø 3.5 mm	140.714	140.715	140.716	140.717	140.718	140.719
ø 3.75 mm	140.720	140.721	140.722	140.723	140.724	140.725
ø 4.0 mm	140.726	140.727	140.728	140.729	140.730	140.731
ø 4.3 mm	140.732	140.733	140.734	140.735	140.736	140.737
ø 5.0 mm	140.738	140.739	140.740	140.741	140.742	140.743

Implant CM Screw Standard

	7.0 mm	8.0 mm	9.0 mm	11.0 mm	13.0 mm	15.0 mm
ø 3.5 mm	140.744	140.745	140.746	140.747	140.748	140.749
ø 3.75 mm	140.750	140.751	140.752	140.753	140.754	140.755
ø 4.0 mm	140.756	140.757	140.758	140.759	140.760	140.761
ø 4.3 mm	140.762	140.763	140.764	140.765	140.766	140.767
ø 5.0 mm	140.768	140.769	-	-	-	-

Implant CM Screw Slim

	7.0 mm	8.0 mm	9.0 mm	11.0 mm	13.0 mm	15.0 mm
ø 3.0 mm	-	-	141.139	141.140	141.141	-

WARNING:

This usage indication should be read before using the Kopp implant system. The Kopp implant system can only be used for its dental purpose, always respecting the dental surgical and work safety laws. CM Screw components should only be used for dental treatment, with the CM Screw Kopp Implant. If the indication and use is not clear, treatment should be suspended until there are clarifications necessary for the proper continuation of the procedure.

Only qualified professionals with knowledge in implantology are authorized to use Kopp products. In addition, the use of products with inappropriate surgical techniques and / or procedures can lead to unsatisfactory results. Therefore, the proper application and use of the products must be respected. Kopp will not assume legal responsibility for product damage resulting from failure to observe the instructions for use.

The following surgical complications were seen after treatment with dental implants:

- Postoperative bleeding.
- Necrosis caused by insufficient cooling.
- Infections.
- Suture dehiscence.
- iatrogenic trauma.
- Implant loss caused by insufficient osseointegration or loosening after installation.
- Periodontal complications caused by insufficient thickness of the mucogingival retention element.
- Prosthetic component stuck or too tight due to excessive torque.
- Fracture of the adaptation screw due to excessive torque.
- Devices in the patient's mouth can be swallowed or aspirated.
- In cases of extremely unfavorable loads (prosthetic overload, high level of bone resorption), rupture of the dental implant body can occur.

PRECAUTIONS

- Always follow the appropriate biosafety requirements before performing any procedure.
- Position the patient in the best possible way to prevent the risk of aspiration of devices as much as possible.
- All devices placed in the patient's mouth must be secure in order to prevent aspiration or swallowing.
- Before the procedure, ensure that all components, instruments and materials are available in the appropriate quantity for the operational execution of the procedure.
- Maximum cutting speeds must be observed beforehand.
- During the bone drilling sequence, maintain the recommended cutting speed, preferably avoiding speed reduction.
- Exchange reusable cutting instruments after 20 consecutive use or when they present any modification to the product that may hinder its functionality.
- All products considered to be for single use should not be reused. Failure to observe this item may cause a risk of infection. Complications resulting from reuse can even cause the loss of the dental implant.

APPLICATIONS: PERIIMPLANT STABILITY

Prosthetic preoperative planning:

In preoperative planning, the length and diameter of the dental implant must be selected to verify that the height and width of the maxillary bone are sufficient for the insertion of the dental implant.

Dental implant installation:

For a successful installation of the dental implant, the dentist must be guided by the knowledge acquired in the field of implant dentistry and make sure that the procedures are followed:

- Carefully carry out the patient's clinical diagnosis.
- Plan the surgery and prepare the patient's medical record. To do this, follow the instructions below:

Dental implant installation:

For a successful installation of the dental implant, the dentist must be guided by the knowledge acquired in the field of implant dentistry and make sure that the procedures are followed:

- Carefully carry out the patient's clinical diagnosis.
- Plan the surgery and prepare the patient's medical record. To do this, follow the instructions below:

In the preoperative period, the patient's treatment must be properly evaluated so that the specifications of the products to be used in the surgical procedure are properly designated by the professional. In addition, the surgical technique must be planned and it is recommended to prepare a controlled protocol of the postoperative period, including the actions taken by the oral surgeon in cases of unexpected interventions during the surgical procedure. In addition, prophylactic treatment should be adopted both preoperatively and in the immediate and late postoperative period as recommended in the classic literature.

Surgical interventions for dental implants should not be considered in circumstances where the patient is not medically healthy. The preoperative clinical evaluation should be performed with caution by the oral surgeon, especially in patients who demonstrate specific concerns or systematic factors that may interfere with the bone or soft tissue healing process.

In addition, the appropriate volume of bone at the surgical site where the dental implant will be implanted should also be analyzed beforehand and during surgery. For a satisfactory treatment, the patient must have an adequate amount of residual bone at the surgical site and the surgeon must select the appropriate height and size of the implant. This clinical evaluation is extremely important for the long-term preservation of the product's functionality. In cases of unfavorable bone volume, bone graft techniques should be considered to be taken beforehand until there is an appropriate amount of bone volume for the installation of the dental implant. Otherwise, in any circumstance where the professional has doubts about the adequate bone volume, the procedure should not be performed.

It is also emphasized that the simultaneous elaboration of a patient's medical record is of utmost importance, having all the information mentioned, as well as the other notes that are relevant for an adequate clinical and surgical follow-up.

- Administer preoperative medication.
- Anesthetize the patient previously from the surgical incision.
- Follow the sequence of suitable drills:

To do this, select the drill sequence according to the height and diameter measurement of the selected dental implant. To clarify the recommended drill sequence, refer to the bone drilling chart available in the Kopp product catalog. When installing implants, do not exceed 800 RPM and always irrigate the surgical area with saline.

- Capture the dental implant according to the instructions available in the catalog, and without contamination of the product in the handling stage. The following guidelines are extremely important:

Dismantle the orange cover tube (primary PVC packaging) in the 180 degree position (open mouth). Luego, insert the Contra Angle adjusted to the Kopp Contra Angle Connection inside the internal hexagon of the dental implant. To do this, apply a light pressure. Also, be careful not to overload your strength. On the contrary, it can cause damage to the product.

Then, keep the orange tape in the same position (180 degrees) and break the gap between the dental implant and the symbol. Finally, manipulate the dental implant with a view to the cliff so that the implant is not damaged and transports the surgical socket.

- Check beforehand the torque recommendations for each product available in the Kopp product catalog. It is extremely important not to exceed the established parameters;
- Install in the direction of the clockwork so that the cervical position of the dental implant reaches the 2 mm infraosseous for Morse tapered interface implants.

Prosthetic sequence:

In the immediate postoperative period, install the implant base, the curator, the prosthetic component according to the chosen loading protocol, respecting its indications and limitations. For more information on prosthetic recommendations, see the Kopp product catalog.

THE DRILLING SEQUENCE RECOMMENDS

	Ø 3.0	Ø 3.5	Ø 3.75	Ø 4.0	Ø 4.3	Ø 5.0
lança	✓	✓	✓	✓	✓	✓
2.0	✓	✓	✓	✓	✓	✓
2.8	✓	✓	✓	✓	✓	✓
3.0		✓	✓	✓	✓	✓
3.3			✓	✓	✓	✓
3.8					✓	✓
4.3						✓

PRODUCT TRACEABILITY:

Traceability label:

This product comes with two adhesive labels that allow its traceability and must be attached to the documents:

- Medical records;
- Documentation required by the patient.

Product traceability requirement:

In any case where the oral surgeon requires further analysis of a particular product, the dentist must provide the patient's history attached with the sticker. On the contrary, the company will be exempted from any legal responsibility in relation to the product. In the case of installing implant-supported prostheses, it is mandatory for the professional to report the technical specifications of the chosen prosthetic component or components.

In addition, the surgeon must take into account beforehand that there is a low percentage of failure in osseointegration due to several factors indicated in the scientific dental literature. Therefore, in none of these occasions the company will not assume responsibility. Furthermore, poor quality or quantity of bone, poor hygiene, infections, and systemic disorders can lead to implant loss in the future.

STERILIZATION:

The dental implant is single use and is supplied sterile, being packaged in a package that offers triple protection. In contrast, not all prosthetic components are sterile.

Important: Sterile products cannot be re-sterilized.

PACKAGING:

Kopp implants are packaged in:

- PVC primary packaging: tube and cap.
- Secondary packaging: blister sealed with surgical grade paper.
- Tertiary packaging - cardboard.

STORAGE CONDITIONS:

Kopp products must be stored, in their original packaging, in a clean and dry place, at room temperature and protected from solar radiation. Keep them out of the reach of children.

VIOLATION OF MATERIALS:

In the event of any damage to the packaging or product, the dental implant must be disposed of immediately and not used under any circumstances.

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REGISTRATION AT THE ANVISA [Federal Government Health Surveillance Agency] Nº: 80204110010

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CNPJ/MF number: 02.967.738/0001-58 State Taxpayer number: 901.793.514-0 -Governmental Authorization for Doing Business: ANVISA/MS

[Federal Government Health Surveillance Agency]: G 2 W 2 8 M X M H 1 3 2 – For more information , please visit: www.koppimplantes.com.br

RESPONSIBLE PROFESSIONAL: GINO KOPP - CRO Nº 7679- PR [Regional Odontology Board- Paraná State].

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